


Product Datasheet

Mouse CCL3 ELISA Kit (orb391042)

Description	Sandwich ELISA for quantitative detection of mouse CCL3 concentrations.
Reactivity	Mouse
Range	7.8 pg/ml - 500 pg/ml
Tested Applications	ELISA
Target	CCL3
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Detection Method: Colorimetric Sample Types: Cell Culture Supernatants, Serum, Plasma
Sample Types	Cell Culture Supernates, Serum, Plasma
Uniprot ID	P10855
Sensitivity	4 pg/ml
Entrez	20302
Expiration Date	Please enquire.



Representative
Standard
Curve