


## Product Datasheet

### Mouse IFN gamma ELISA Kit (orb391083)

<b>Description</b>	Sandwich ELISA for quantitative detection of mouse IFN gamma concentrations.
<b>Reactivity</b>	Mouse
<b>Range</b>	7.8 pg/ml - 500 pg/ml
<b>Tested Applications</b>	ELISA
<b>Target</b>	IFNG
<b>Storage</b>	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
<b>Note</b>	For research use only
<b>Application notes</b>	Detection Method: Colorimetric Sample Types: Cell Culture Supernatants, Serum, Plasma, Tissue Homogenates
<b>Sample Types</b>	Cell Culture Supernates, Serum, Plasma, Tissue Homogenates
<b>Uniprot ID</b>	<b>P01580</b>
<b>Sensitivity</b>	4 pg/ml
<b>Entrez</b>	<b>15978</b>
<b>Expiration Date</b>	Please enquire.



Representative  
Standard  
Curve