


Product Datasheet

Mouse IGFBP2 ELISA Kit (orb391088)

Description	Sandwich ELISA for quantitative detection of mouse IGFBP2 concentrations.
Reactivity	Mouse
Range	312 pg/ml - 20000 pg/ml
Tested Applications	ELISA
Target	IGFBP2
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Detection Method: Colorimetric Sample Types: Cell Culture Supernatants, Serum, Plasma
Sample Types	Cell Culture Supernates, Serum, Plasma
Uniprot ID	P47877
Sensitivity	120 pg/ml
Entrez	16008
Expiration Date	Please enquire.



Representative
Standard
Curve