


Product Datasheet

Rat IFN gamma ELISA Kit (orb391157)

Description	Sandwich ELISA for quantitative detection of rat IFN gamma concentrations.
Reactivity	Rat
Range	31.2 pg/ml - 2000 pg/ml
Tested Applications	ELISA
Target	IFNG
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Detection Method: Colorimetric Sample Types: Cell Culture Supernatants, Serum, Plasma, Tissue Homogenates
Sample Types	Cell Culture Supernates, Serum, Plasma, Tissue Homogenates
Uniprot ID	P01581
Sensitivity	7 pg/ml
Entrez	25712
Expiration Date	Please enquire.



Representative
Standard
Curve