


Product Datasheet

Rat IGF2 ELISA Kit (orb391159)

Description	Sandwich ELISA for quantitative detection of rat IGF2 concentrations.
Reactivity	Rat
Range	46.9 pg/ml - 3000 pg/ml
Tested Applications	ELISA
Target	IGF2
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Detection Method: Colorimetric Sample Types: Cell Culture Supernatants, Serum, Plasma
Sample Types	Cell Culture Supernates, Serum, Plasma
Uniprot ID	P01346
Sensitivity	25 pg/ml
Entrez	24483
Expiration Date	Please enquire.



Representative Standard Curve