
Product Datasheet

Mouse ADIPOQ protein (orb391203)

Description

Recombinant mouse ADIPOQ protein is

and the target gene en...

Reactivity

Mouse

Endotoxins

Endotoxin: Less than 0.1 ng/μg (1 IEU/μg) as determined by LAL test.

Tested Applications

ELISA, MS, SDS-PAGE, WB

Target

ADIPOQ

Preservatives

Lyophilized from a 0.2 μM filtered solution of PBS, pH 7.4.

Form/Appearance

Lyophilized from a 0.2 μM filtered solution of PBS, pH 7.4.

Storage

Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 2-8°C for 2-7 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.

Note

For research use only

Application notes

Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100 μg/ml. Dissolve the lyophilized protein in 1X PBS. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

Protein Sequence

MGSSHHHHHH SSGLVPRGSH MEDDVTTTEE
 LAPALVPPPK GTCAGWMAGI PGHPGHNGTP
 GRDGRDGTPG EKGEKGDAGL LGPKGETGDV
 GMTGAEGPRG FPGTPGRKGE PGEAAYVYRS
 AFSVGLETRV TVPNVPIRFT KIFYNQNHY
 DGSTGKFYCN IPGLYYFSYH ITVYMKDKVK
 SLFKKDKAVL FTYDQYQEKV VDQASGSVLL
 HLEVGDQVWL QVYGDGDHNG LYADNVNDST
 FTGFLLYHDT N

Purity

Greater than 95% as determined by reducing SDS-PAGE.

Source

E. coli

Uniprot ID
[Q60994](#)
Entrez
[11450](#)