

Product Datasheet

Mouse ADIPOQ protein (orb391204)

Description

Recombinant mouse ADIPOQ protein is

system and the target gene ...

Reactivity

Mouse

Endotoxins

Endotoxin: Less than 0.1 ng/μg (1 IEU/μg) as determined by LAL test.

Tested Applications

ELISA, MS, SDS-PAGE, WB

Target

ADIPOQ

Preservatives

Lyophilized from a 0.2 μM filtered solution of 20mM PB, 150mM NaCl, pH 7.4.

Form/Appearance

Lyophilized from a 0.2 μM filtered solution of 20mM PB, 150mM NaCl, pH 7.4.

Storage

Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 2-8°C for 2-7 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.

Note

For research use only

Application notes

Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100 μg/ml. Dissolve the lyophilized protein in 1X PBS. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

Protein Sequence

EDDVTTTEEL APALVPPPKG TCAGWMAGIP
 GHPGHNGTPG RDGRDGTPE KGEKGDAGLL
 GPKGETGDVG MTGAEGPRGF PGTPGRKGEP
 GEAAVYRSA FSVGLETRVT VPNVPIRFTK
 IFYNQQNHVD GSTGKFCNI PGLYYFSYHI
 TVYMKDVKVS LFKKDKAVLF TYDQYQEKNV
 DQASGSVLLH LEVGDQVWLQ VYGDGDHNGL
 YADNVNDSTF TGFLLYHDTN HHHHHH

Purity

Greater than 95% as determined by reducing SDS-PAGE.

Source

Human cells

Uniprot ID
[Q60994](#)
Entrez
[11450](#)
Expiration Date

6 months from date of receipt.