

---

## Product Datasheet

### Human C19ORF10 protein (orb391320)

**Description**

Recombinant human C19ORF10 protein is

and the target gene ...

**Reactivity**

Human

**Endotoxins**

Endotoxin: Less than 0.1 ng/μg (1 IEU/μg) as determined by LAL test.

**Tested Applications**

ELISA, MS, SDS-PAGE, WB

**Target**

C19ORF10

**Preservatives**

Lyophilized from a 0.2 μM filtered solution of 20mM PB, 150mM NaCl, pH 7.4.

**Form/Appearance**

Lyophilized from a 0.2 μM filtered solution of 20mM PB, 150mM NaCl, pH 7.4.

**Storage**

Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 2-8°C for 2-7 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.

**Note**

For research use only

**Application notes**

Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100 μg/ml. Dissolve the lyophilized protein in 1X PBS. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

**Protein Sequence**

MGSSHHHHHH SSGLVPRGSH MSEPTTVAFD  
 VRPGGVVHSF SHNVGPGDKY TCMFTYASQG  
 GTNEQWQMSL GTSEDHQHFT CTIWRPQGKS  
 YLYFTQFKAE VRGAEIEYAM AYSKAAFERE  
 SDVPLKTEEF EVTKTAVAGR PGAFKAELSK  
 LVIVAKASRT EL

**Purity**

Greater than 95% as determined by reducing SDS-PAGE.

**Source**

E. coli

**Uniprot ID**
[Q969H8](#)
**Entrez**
[56005](#)
**Expiration Date**

6 months from date of receipt.