

Product Datasheet

Human CCL27 protein (orb391370)

Description

Recombinant human CCL27 protein is

glycosylated polypept...

Reactivity

Human

Endotoxins

Endotoxin: Less than 0.1 ng/μg (1 IEU/μg) as determined by LAL test.

Tested Applications

ELISA, MS, SDS-PAGE, WB

Target

CCL27

Preservatives

Lyophilized from a 0.2 μM filtered solution of 20mM PB, 500mM NaCl, pH 7.5

Form/Appearance

Lyophilized from a 0.2 μM filtered solution of 20mM PB, 500mM NaCl, pH 7.5

Storage

Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 2-8°C for 2-7 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.

Note

For research use only

Application notes

Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100 μg/ml. Dissolve the lyophilized protein in 1X PBS. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

Protein Sequence

 FLLPPSTACC TQLYRKPLSD KLLRKVIQVE
 LQEADGDCHL QAFVLHLAQR SICIHPQNPS
 LSQWFEHQER KLHGTLPKLN FGMLRKMG

Purity

Greater than 95% as determined by SEC-HPLC and reducing SDS-PAGE.

Source

E. coli

Uniprot ID
[Q9Y4X3](#)
Entrez
[10850](#)
Expiration Date

6 months from date of receipt.