

---

## Product Datasheet

### Human CD300C protein (orb391403)

**Description**

Recombinant human CD300C protein is

system and the target gene ...

**Reactivity**

Human

**Endotoxins**

Endotoxin: Less than 0.1 ng/μg (1 IEU/μg) as determined by LAL test.

**Tested Applications**

ELISA, MS, SDS-PAGE, WB

**Target**

CD300C

**Preservatives**

Lyophilized from a 0.2 μM filtered solution of 20mM PB, 150mM NaCl, pH 7.4.

**Form/Appearance**

Lyophilized from a 0.2 μM filtered solution of 20mM PB, 150mM NaCl, pH 7.4.

**Storage**

Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 2-8°C for 2-7 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.

**Note**

For research use only

**Application notes**

Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100 μg/ml. Dissolve the lyophilized protein in 1X PBS. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

**Protein Sequence**

GYFPLSHPMT VAGPVGGSLS VQCRYEKEHR  
 TLNKFWRPP QILRCDKIVE TKGSAGKRNG  
 RVSIRDSPAN LSFTVTLENL TEEDAGTYWC  
 GVDTPWLRDF HDPIVEVEVS VFPAGTTTAS  
 SPQSSMGTS G PPTKLPVHTW PSVTRKDSPE  
 PSPHPGSLFS NVRVDDIEGR MDEPKSCDKT  
 HTCPPCPAPE LLGGPSVFLF PPKPKDTLMI  
 SRTPEVTCVV VDVSHEDPEV KFNWYVDGVE  
 VHNAKTKPRE EQYNSTYRVV SVLTVLHQDW  
 LNGKEYKCKV SNKALPAPIE KTISKAKGQP  
 REPQVYTLPP SREEMTKNQV SLTCLVKGFY  
 PSDIAVEWES NGQPENNYKT TPPVLDSGDS  
 FFLYSKLTVD KSRWQQGNVF SCSVMHEALH  
 NHYTQKLSL SPGK

**Purity**

Greater than 95% as determined by reducing SDS-PAGE.