

---

## Product Datasheet

### Human CLIC3 protein (orb391490)

**Description**

Recombinant human CLIC3 protein is

and the target gene enc...

**Reactivity**

Human

**Endotoxins**

Endotoxin: Less than 0.1 ng/μg (1 IEU/μg) as determined by LAL test.

**Tested Applications**

ELISA, MS, SDS-PAGE, WB

**Target**

CLIC3

**Preservatives**

Liquid in a 0.2 μM filtered solution of 10mM Tris, 0.1%Triton100, pH 8.0.

**Form/Appearance**

Liquid in a 0.2 μM filtered solution of 10mM Tris, 0.1%Triton100, pH 8.0.

**Storage**

Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

**Note**

For research use only

**Application notes**

Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100 μg/ml. Dissolve the lyophilized protein in 1X PBS. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

**Protein Sequence**

 MAETKLQLFV KASEDGESVG HCPSCQRLFM  
 VLLLKGV PFT LTTVDTRRSP DVLKDFAPGS  
 QLPILLYDSD AKTDTLQIED FLEETLGPPD  
 FPSLAPRYRE SNTAGNDVFH KFSAFIKNPV  
 PAQDEALYQQ LLRALARLDS YLRAPLEHEL  
 AGEPQLRESR RRFLDGDRLT LADCSLLPKL  
 HIVDTVCAHF RQAPIPAELR GVERRYLDSAM  
 QEKEFKYTCP HSAEILAAAYR PAVHPRLEHH  
 HHHH

**Purity**

Greater than 95% as determined by reducing SDS-PAGE.

**Source**

E. coli

**Uniprot ID**
**095833**
**Entrez**
**9022**
**Expiration Date**

6 months from date of receipt.