
Product Datasheet

Human FLRT1 protein (orb391750)

Description

Recombinant human FLRT1 protein is

system and the target gene e...

Reactivity

Human

Endotoxins

Endotoxin: Less than 0.1 ng/μg (1 IEU/μg) as determined by LAL test.

Tested Applications

ELISA, MS, SDS-PAGE, WB

Target

FLRT1

Preservatives

Lyophilized from a 0.2 μM filtered solution of 20mM PB, 150mM NaCl, pH 7.2.

Form/Appearance

Lyophilized from a 0.2 μM filtered solution of 20mM PB, 150mM NaCl, pH 7.2.

Storage

Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 2-8°C for 2-7 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.

Note

For research use only

Application notes

Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100 μg/ml. Dissolve the lyophilized protein in 1X PBS. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

Protein Sequence

```

IDSTTCPSVC RCDNGFIYCN DRGLTSIPAD
IPDDATTLYL QNNQINNAGI PQDLKTKVNV
QVIYLYENDL DEFPINLPRS LRELHLQDNN
VRTIARDSLA RIPPLEKLHL DDNSVSTVSI
EEDAFADSKQ LKLLFLSRNH LSSIPGLPH
TLEELRLDDN RISTIPLHAF KGLNSLRRLV
LDGNLLANQR IADDTFSRLQ NLTELSLVRN
SLAAPPLNLP SAHLQKLYLQ DNAISHIPYN
TLAKMRELER LDLSNNNLTT LPRGLFDDLG
NLAQLLLRRN PWFCGCNLMW LRDWVKARAA
VVNVRGLMCQ GPEKVRGMAI KDITSEMDEC
FETGPQGGVA NAAAKTTASN HASATTPQGS
LFTLKAKRPG LRLPDSNIDY PMATGDGAKT
LAIHVKALTA DSIRITWKAT LPASSFRLSW
LRLGHSPAVG SITETLVQGD KTEYLLTALE
PKSTYIICMV TMETSNAVVA DETPVCARAE
TADSYGPTTT LNQEQNAGPM ASLPVDHHHH
HH
    
```

Biorbyt Ltd.

7 Signet Court, Swann's Road, Cambridge, CB5 8LA, United Kingdom

 Email: info@biorbyt.com | Phone: +44 (0) 1223 859-353 | Fax: +44 (0)1223 280 240

Biorbyt LLC.

68 TW Alexander Drive
Research Triangle Park
Durham, North Carolina
27709, United States

 Email: info@biorbyt.com | Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558