
Product Datasheet

Human GFRA2 protein (orb391797)

Description

Recombinant human GFRA2 protein is

system and the target gene e...

Reactivity

Human

Endotoxins

Endotoxin: Less than 0.1 ng/μg (1 IEU/μg) as determined by LAL test.

Tested Applications

ELISA, MS, SDS-PAGE, WB

Target

GFRA2

Preservatives

Lyophilized from a 0.2 μM filtered solution of 20mM Tris, 150mM NaCl, pH 8.0.

Form/Appearance

Lyophilized from a 0.2 μM filtered solution of 20mM Tris, 150mM NaCl, pH 8.0.

Storage

Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 2-8°C for 2-7 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.

Note

For research use only

Application notes

Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100 μg/ml. Dissolve the lyophilized protein in 1X PBS. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

Protein Sequence

SPSSLQGPEL HGWRPPVDCV RANELCAAES
 NCSSRYRTL R QCLAGRDRNT MLANKECQAA
 LEVLQESPLY DCRCKRGMKK ELQCLQIYWS
 IHLGLTEGEE FYEASPYEPV TSRLSDIFRL
 ASIFSGTGAD PVVSAKSNHC LDAKACNLN
 DNCKKLRSSY ISICNREISP TERCNRRKCH
 KALRQFFDRV PSEYTYRMLF CSCQDQACAE
 RRRQTILPSC SYEDKEKPNC LDLRGVCRTD
 HLCRSRLADF HANCRASYQT VTSCPADNYQ
 ACLGSYAGMI GFDMPNYVD SSPTGIVVSP
 WCSCRGSGNM EEECEKFLRD FTENPCLRNA
 IQAFGNGTDV NVSPKGPSFQ ATQAPRVEKT
 PSLPDDLSDS TSLGTSVITT CTSVQEQGLK
 ANNSKELSMC FTELTTNIIP GSNKVIKPN
 VDDIEGRMDE PKSCDKTHTC PPCPAPELLG
 GPSVFLFPPK PKDTLMISRT PEVTCVVVDV
 SHEDPEVKFN WYVDGVEVHN AKTKPREEQY
 NSTYRVVSVI TVI HODWING KEYKCKVSNK

Biorbyt Ltd.

7 Signet Court, Swann's Road, Cambridge, CB5 8LA, United Kingdom

 Email: info@biorbyt.com | Phone: +44 (0) 1223 859-353 | Fax: +44 (0)1223 280 240

Biorbyt LLC.

68 TW Alexander Drive
Research Triangle Park
Durham, North Carolina
27709, United States

 Email: info@biorbyt.com | Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558

Purity

Greater than 95% as determined by reducing