
Product Datasheet

Human IL10RB protein (orb391911)

Description

Recombinant human IL10RB protein is

system and the target gene ...

Reactivity

Human

Endotoxins

Endotoxin: Less than 0.1 ng/μg (1 IEU/μg) as determined by LAL test.

Tested Applications

ELISA, MS, SDS-PAGE, WB

Target

IL10RB

Preservatives

Lyophilized from a 0.2 μM filtered solution of PBS, pH 7.4.

Form/Appearance

Lyophilized from a 0.2 μM filtered solution of PBS, pH 7.4.

Storage

Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 2-8°C for 2-7 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.

Note

For research use only

Application notes

Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100 μg/ml. Dissolve the lyophilized protein in 1X PBS. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

Protein Sequence

MVPPPENVRM NSVNFKNILQ WESPAFAKGN
 LTFTAQYLSY RIFQDKCMNT TLTECDFSSL
 SKYGDHTLRV RAEFADEHSD WVNITFCPVD
 DTIIGPPGMQ VEVLADSLHM RFLAPKIENE
 YETWTMKNVY NSWTYNVQYW KNGTDEKFQI
 TPQYDFEVLN NLEPWTTYCV QVRGFLPDRN
 KAGEWSEPVC EQTTHDETVP SVDDIEGRMD
 EPKSCDKTHT CPPCPAPELL GGPSVFLFPP
 KPKDTLMISR TPEVTCVVVD VSHEDPEVKF
 NWYVDGVEVH NAKTKPREEQ YNSTYRVVSV
 LTVLHQDWLN GKEYKCKVSN KALPAIEKT
 ISKAKGQPRE PQVYTLPPSR EEMTKNQVSL
 TCLVKGFYPS DIAVEWESNG QPENNYKTPP
 PVLDSGDGFF LYSKLTVDKS RWQQGNVFSC
 SVMHEALHNN YTQKLSLSLSP GK

Purity

Greater than 95% as determined by reducing SDS-PAGE.

Biorbyt Ltd.

7 Signet Court, Swann's Road, Cambridge, CB5 8LA, United Kingdom

 Email: info@biorbyt.com | Phone: +44 (0) 1223 859-353 | Fax: +44 (0)1223 280 240

Biorbyt LLC.

68 TW Alexander Drive
Research Triangle Park
Durham, North Carolina
27709, United States

 Email: info@biorbyt.com | Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558