
Product Datasheet

Human IL17RA protein (orb391935)

Description

Recombinant human IL17RA protein is

system and the target gene ...

Reactivity

Human

Endotoxins

Endotoxin: Less than 0.1 ng/μg (1 IEU/μg) as determined by LAL test.

Tested Applications

ELISA, MS, SDS-PAGE, WB

Target

IL17RA

Preservatives

Lyophilized from a 0.2 μM filtered solution of PBS, pH 7.4.

Form/Appearance

Lyophilized from a 0.2 μM filtered solution of PBS, pH 7.4.

Storage

Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 2-8°C for 2-7 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.

Note

For research use only

Application notes

Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100 μg/ml. Dissolve the lyophilized protein in 1X PBS. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

Protein Sequence

LRLLDHRALV CSQPGLNCTV KNSTCLDDSW
 IHPRNLTSS PKDLQIQLHF AHTQQGDLFP
 VAHIEWTLQT DASILYLEGA ELSVLQNTN
 ERLCVRFEFL SKLRHHHRRW RFTFSHFVVD
 PDQEYEVTVH HLPKPIPDGD PNHQSKNFLV
 PDCEHARMKV TTPCMSSGSL WDPNITVETL
 EAHQLRVSFT LWNETHYQI LLTSFPHMEN
 HSCFEHMHHI PAPRPEEFHQ RSNVTLTLRN
 LKGCCRHQVQ IQPFSSCLN DCLRHSATVS
 CPEMPDTPEP IPDYMLWVD DIEGRMDEPK
 SCDKTHTCPP CPAPELLGGP SVFLFPPKPK
 DTLMISRTPE VTCVVVDVSH EDPEVKFNWY
 VDGVEVHNAK TKPREEQYNS TYRVVSVLTV
 LHQDWLNGKE YKCKVSNKAL PAPIEKTISK
 AKGQPREPQV YTLPPSREEM TKNQVSLTCL
 VKGFYPSDIA VEWESNGQPE NNYKTTTPVL
 DSDGSFFLYS KLTVDKSRWQ QGNVFSQSV
 HEALHNHYTO KSI SI SPGK