
Product Datasheet

Mouse IL1 alpha protein (orb391944)

Description

Recombinant mouse IL1 alpha protein is

glycosylated poly...

Reactivity

Mouse

Endotoxins

Endotoxin: Less than 0.1 ng/μg (1 IEU/μg) as determined by LAL test.

Tested Applications

ELISA, MS, SDS-PAGE, WB

Target

IL1A

Preservatives

Lyophilized from a 0.2 μM filtered solution of 50mM Tris-HCl, 200mM NaCl, pH 8.0.

Form/Appearance

Lyophilized from a 0.2 μM filtered solution of 50mM Tris-HCl, 200mM NaCl, pH 8.0.

Storage

Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 2-8°C for 2-7 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.

Note

For research use only

Application notes

Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100 μg/ml. Dissolve the lyophilized protein in 1X PBS. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

Protein Sequence

 SAPYTYQSDL RYKLMKLVQR KFMNDSLNDQ
 TIYQDVVKHY LSTTWLNDLQ QEVKFDYAY
 SSGGDDSKYP VTLKISDSQL FVSAQGEDQP
 VLLKELPETP KLITGSETDL IFFWKSINSK
 NYFTSAAYPE LFIATKEQSR VHLARGLPSM
 TDFQIS

Purity

Greater than 95% as determined by SEC-HPLC and reducing SDS-PAGE.

Source

E. coli

Uniprot ID
[P01582](#)
Entrez
[16175](#)
Expiration Date

6 months from date of receipt.