
Product Datasheet

Human IL1RA protein (orb391952)

Description

Recombinant human IL1RA protein is

and the target gene enc...

Reactivity

Human

Endotoxins

Endotoxin: Less than 0.1 ng/μg (1 IEU/μg) as determined by LAL test.

Tested Applications

ELISA, MS, SDS-PAGE, WB

Target

IL1RN

Preservatives

Lyophilized from a 0.2 μM filtered solution of 50mM Tris-HCl, 200mM NaCl, pH 7.5.

Form/Appearance

Lyophilized from a 0.2 μM filtered solution of 50mM Tris-HCl, 200mM NaCl, pH 7.5.

Storage

Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 2-8°C for 2-7 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.

Note

For research use only

Application notes

Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100 μg/ml. Dissolve the lyophilized protein in 1X PBS. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

Protein Sequence

MRPSGRKSSK MQAFRIWDVN QKTFYLRNNQ
 LVAGYLQGN VNLEEKIDVV PIEPHALFLG
 IHGGKMCLSC VKSGDETRLQ LEAVNITDLS
 ENRKQDKRFA FIRSDSGPTT SFESAACPGW
 FLCTAMEADQ PVSLTNMPDE GVMVTKFYFQ
 EDE

Purity

Greater than 95% as determined by RP-HPLC and reducing SDS-PAGE.

Source

E. coli

Uniprot ID
[P18510](#)
Entrez
[3557](#)
Expiration Date

6 months from date of receipt.