

---

## Product Datasheet

### Human IL6 protein (orb391992)

**Description**

Recombinant human IL6 protein is produced

gene encodi...

**Reactivity**

Human

**Endotoxins**

Endotoxin: Less than 0.1 ng/μg (1 IEU/μg) as determined by LAL test.

**Tested Applications**

ELISA, MS, SDS-PAGE, WB

**Target**

IL6

**Preservatives**

Lyophilized from a 0.2 μM filtered solution of 20mM HAc-NaAc, 150mM NaCl, pH 5.0

**Form/Appearance**

Lyophilized from a 0.2 μM filtered solution of 20mM HAc-NaAc, 150mM NaCl, pH 5.0

**Storage**

Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 2-8°C for 2-7 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.

**Note**

For research use only

**Application notes**

Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100 μg/ml. Dissolve the lyophilized protein in 1X PBS. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

**Protein Sequence**

PVPPGEDSKD VAAPHRQPLT SSERIDKQIR  
 YILDGISALR KETCNKSNMC ESSKEALAEN  
 NLNLPKMAEK DGCFSQGFNE ETCLVKIITG  
 LLEFEVYLEY LQNRFESSEE QARAVQMSTK  
 VLIQFLQKKA KNLDAITTPD PTTNASLLTK  
 LQAQNQWLQD MTTHLILRSF KEFLQSSLRA  
 LRQM

**Purity**

Greater than 98% as determined by RP-HPLC and reducing SDS-PAGE.

**Source**

E. coli

**Uniprot ID**
[P05231](#)
**Entrez**
[3569](#)
**Expiration Date**

6 months from date of receipt.