
Product Datasheet

Human PACSIN2 protein (orb392230)

Description

Recombinant human PACSIN2 protein is

system and the target gene...

Reactivity

Human

Endotoxins

Endotoxin: Less than 0.1 ng/μg (1 IEU/μg) as determined by LAL test.

Tested Applications

ELISA, MS, SDS-PAGE, WB

Target

PACSIN2

Preservatives

Lyophilized from a 0.2 μM filtered solution of 20mM PB, 150mM NaCl, pH 7.4.

Form/Appearance

Lyophilized from a 0.2 μM filtered solution of 20mM PB, 150mM NaCl, pH 7.4.

Storage

Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 2-8°C for 2-7 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.

Note

For research use only

Application notes

Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100 μg/ml. Dissolve the lyophilized protein in 1X PBS. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

Protein Sequence

```

MSVTYDDSVG VEVSSDSFWE VGNYKRTVKR
IDDGHRLCD LMNCLHERAR IEKAYAQLT
EWARRWRQLV EKG PQYGTVE KAWMAFMSEA
ERVSELHLEV KASLMNDDFE KIKNWQKEAF
HKQMMGGFKE TKEAEDGFRK AQPWAKKLL
EVEAAKKAHH AACKEEKLA I SREANSKADP
SLNPEQLKKL QDKIEKCKQD VLKTKKEYEK
SLKELDQGTP QYMNMEQVF EQCQQFEEKR
LRFFREV LLE VQKHLDSL NV AGYKAIYHDL
EQSIRAADAV EDLRWFRANH GPGMAMNWPQ
FEEWSADLNR TLRREKKKA TDGVTLTGIN
QTGDQSLPSK PSSTLNVP SN PAQSAQSQSS
YNPFEDEDDT GSTVSEKDDT KAKNVSSYEK
TQSYPTDWS D DESNNPFSST DANGDSNPF D
DDATSGTEVR VRALYD YEGQ EHD ELSFKAG
DELTKMEDED EQGWCKGR LD NGQVGLYPAN
YVEAIQVDHH HHHH
  
```