

Product Datasheet

Human PAFAH1B2 protein (orb392232)

Description

Recombinant human PAFAH1B2 protein is

and the target gene e...

Reactivity

Human

Endotoxins

Endotoxin: Less than 0.1 ng/μg (1 IEU/μg) as determined by LAL test.

Tested

ELISA, MS, SDS-PAGE, WB

Applications
Target

PAFAH1B2

Preservatives

Liquid in a 0.2 μM filtered solution of 20mM Tris, 150mM NaCl, pH 8.0.

Form/Appearance

Liquid in a 0.2 μM filtered solution of 20mM Tris, 150mM NaCl, pH 8.0.

Storage

Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Note

For research use only

Application notes

Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100 μg/ml. Dissolve the lyophilized protein in 1X PBS. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

Protein Sequence

MSQGDSPAA IPHAAEDIQG DDRWMSQHNR
FVLDCCKDEP DVLFVGDMSV QLMQQYEIWR
ELFSPLHALN FGIGGDTTRH VLWRLKNGEL
ENIKPKVIVV WVGTTNNHENT AEEVAGGIEA
IVQLINTRQP QAKIIVLGLL PRGEKPNPLR
QKNAKVNQLL KVSPLKLANV QLLDTDGGFV
HSDGAISCHD MFDLHLTGG GYAKICKPLH
ELIMQLEET PEEKQTTIAV EHHHHHHH

Purity

Greater than 95% as determined by reducing SDS-PAGE.

Source

E. coli

Uniprot ID

P68402

Entrez

5049

Expiration Date

6 months from date of receipt.