
Product Datasheet

Human RSP01 protein (orb392401)

Description

Recombinant human RSPO1 protein is

system and the target gene e...

Reactivity

Human

Endotoxins

Endotoxin: Less than 0.1 ng/μg (1 IEU/μg) as determined by LAL test.

Tested Applications

ELISA, MS, SDS-PAGE, WB

Target

RSPO1

Preservatives

Lyophilized from a 0.2 μM filtered solution of PBS, pH 7.4.

Form/Appearance

Lyophilized from a 0.2 μM filtered solution of PBS, pH 7.4.

Storage

Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 2-8°C for 2-7 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.

Note

For research use only

Application notes

Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100 μg/ml. Dissolve the lyophilized protein in 1X PBS. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

Protein Sequence

 RISAEGSQAC AKGCELCSEV NGCLKCSPKL
 FILLERNDIR QVGVCLPSCP PGYFDARNPD
 MNKCIKCKIE HCEACFSHNF CTKCKEGLYL
 HKGRCYPACP EGSSAANGTM ECSSPAQCEM
 SEWSPWGPCS KKQQLCGFRR GSEERTRRVL
 HAPVGDHAAC SDTKETRRCT VRRVPCPEGQ
 KRRKGGQGRR ENANRNLARK ESKEAGAGSR
 RRRKGQQQQQ QGTVGPLTSA GPAVDHHHHH
 H

Purity

Greater than 95% as determined by reducing SDS-PAGE.

Source

Human cells

Uniprot ID
[Q2MKA7](#)
Entrez
[284654](#)