

---

## Product Datasheet

### Human SUMO1 protein (orb392539)

**Description**

Recombinant human SUMO1 protein is

and the target gene enc...

**Reactivity**

Human

**Endotoxins**

Endotoxin: Less than 0.1 ng/μg (1 IEU/μg) as determined by LAL test.

**Tested Applications**

ELISA, MS, SDS-PAGE, WB

**Target**

SUMO1

**Preservatives**

Lyophilized from a 0.2 μM filtered solution of 50mM Tris-HCl, 100mM NaCl, 1mM DTT, pH 8.5.

**Form/Appearance**

Lyophilized from a 0.2 μM filtered solution of 50mM Tris-HCl, 100mM NaCl, 1mM DTT, pH 8.5.

**Storage**

Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 2-8°C for 2-7 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.

**Note**

For research use only

**Application notes**

Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100 μg/ml. Dissolve the lyophilized protein in 1X PBS. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

**Protein Sequence**

MGSSHHHHHH SSGLVPRGSH MSDQEAKPST  
EDLGDKKEGE YIKLKVIGQD SSEIHFVKVM  
TTHLKKLKES YCQRQGVPMN SLRFLFEGQR  
IADNNTPKEL GMEEEDVIEV YQEQTGGHST V

**Purity**

Greater than 95% as determined by SEC-HPLC and reducing SDS-PAGE.

**Source**

E. coli

**Uniprot ID**

[P63165](#)

**Entrez**

[7341](#)

**Expiration Date**

6 months from date of receipt.