
Product Datasheet

Human TGF beta Receptor 2 protein (orb392558)

Description

Recombinant human TGF beta Receptor

expression system and the...

Reactivity

Human

Endotoxins

Endotoxin: Less than 0.1 ng/μg (1 IEU/μg) as determined by LAL test.

Tested Applications

ELISA, MS, SDS-PAGE, WB

Target

TGFBR2

Preservatives

Lyophilized from a 0.2 μM filtered solution of 20mM PB, 150mM NaCl, pH 7.2.

Form/Appearance

Lyophilized from a 0.2 μM filtered solution of 20mM PB, 150mM NaCl, pH 7.2.

Storage

Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 2-8°C for 2-7 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.

Note

For research use only

Application notes

Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100 μg/ml. Dissolve the lyophilized protein in 1X PBS. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

Protein Sequence

TIPPHVQKSV NNDMIVTDNN
 GAVKFPQLCK FCDVRFSTCD
 NQKSCMSNCS ITSICEKPQE
 VCVAVWRKND ENITLETVCH DPKLPYHDFI
 LEDAASPKCI MKEKKKPGET FFMCS CSSDE
 CNDNIIFSEE YNTSNPDVDD IEGRMDEPKS
 CDKTHTCPPC PAPELLGGPS VFLFPPKPKD
 TLMISRTPEV TCVVVDVSHE
 DPEVKFNWYV DGVEVHNAKT
 KPREEQYNST YRVVSVLTVL
 HQDWLNGKEY KCKVSNKALP APIEKTISKA
 KGQPREPQVY TLPPSREEMT
 KNQVSLTCLV KGFYPSDIAV
 EWESNGQPEN NYKTTTPVLD
 SDGSEFI YSK I TVDKSRWOO

Biorbyt Ltd.

7 Signet Court, Swann's Road, Cambridge, CB5 8LA, United Kingdom

 Email: info@biorbyt.com | Phone: +44 (0) 1223 859-353 | Fax: +44 (0)1223 280 240

Biorbyt LLC.

68 TW Alexander Drive
Research Triangle Park
Durham, North Carolina
27709, United States

 Email: info@biorbyt.com | Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558