
Product Datasheet

Human TNFRSF10C protein (orb392586)

Description

Recombinant human TNFRSF10C protein is

system and the target ge...

Reactivity

Human

Endotoxins

Endotoxin: Less than 0.1 ng/μg (1 IEU/μg) as determined by LAL test.

Tested Applications

ELISA, MS, SDS-PAGE, WB

Target

TNFRSF10C

Preservatives

Lyophilized from a 0.2 μM filtered solution of 20mM PB, 150mM NaCl, pH 7.2.

Form/Appearance

Lyophilized from a 0.2 μM filtered solution of 20mM PB, 150mM NaCl, pH 7.2.

Storage

Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 2-8°C for 2-7 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.

Note

For research use only

Application notes

Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100 μg/ml. Dissolve the lyophilized protein in 1X PBS. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

Protein Sequence

 ATTARQEEVP QQTVPQQQR HSFKGEECPA
 GSHRSEHTGA CNPCTEGVDY TNASNNEPSC
 FPCTVCKSDQ KHKSSCTMTR DTVCQCKEGT
 FRNENSPEMC RKCSRCPGSGE VQVSNCTSWD
 DIQCVVEEFGA NATVETPAAE ETMNTSPGTP
 APAAEETMNT SPGTPAPAAE ETMTTSPGTP
 APAAEETMTT SPGTPAVDHH HHHH

Purity

Greater than 95% as determined by SEC-HPLC and reducing SDS-PAGE.

Source

Human cells

Uniprot ID
[O14798](#)
Entrez
[8794](#)
Expiration Date

6 months from date of receipt.