

---

## Product Datasheet

### Human TNFRSF1B protein (orb392598)

**Description**

Recombinant human TNFRSF1B protein is

system and the target gen...

**Reactivity**

Human

**Endotoxins**

Endotoxin: Less than 0.1 ng/μg (1 IEU/μg) as determined by LAL test.

**Tested Applications**

ELISA, MS, SDS-PAGE, WB

**Target**

TNFRSF1B

**Preservatives**

Lyophilized from a 0.2 μM filtered solution of PBS, pH 7.4.

**Form/Appearance**

Lyophilized from a 0.2 μM filtered solution of PBS, pH 7.4.

**Storage**

Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 2-8°C for 2-7 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.

**Note**

For research use only

**Application notes**

Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100 μg/ml. Dissolve the lyophilized protein in 1X PBS. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

**Protein Sequence**

 KKKPLCLQRE AKVPHLPADK ARGTQGPEQQ  
 HLLITAPSSS SSSLESSASA LDRRAPTRNQ  
 PQAPGVEASG AGEARASTGS SDSSPGGHGT  
 QVNVTCIVNV CSSSDHSSQC SSQASSTMGD  
 TDSSPSESPK DEQVPFSKEE CAFRSQLETP  
 ETLLGSTEEL PLPLGVDPAG MKPSVDHHHH  
 HH

**Purity**

Greater than 95% as determined by reducing SDS-PAGE.

**Source**

Human cells

**Uniprot ID**
**P20333**
**Entrez**
**7133**
**Expiration Date**

6 months from date of receipt.