

Product Datasheet

Human UBE2L6 protein (orb392665)

Description

Recombinant human UBE2L6 protein is

and the target gene en...

Reactivity

Human

Endotoxins

Endotoxin: Less than 0.1 ng/μg (1 IEU/μg) as determined by LAL test.

Tested

Applications

ELISA, MS, SDS-PAGE, WB

Target

UBE2L6

Preservatives

Liquid in a 0.2 μM filtered solution of 50mM HEPES, 100mM NaCl, pH 8.0.

Form/Appearance

Liquid in a 0.2 μM filtered solution of 50mM HEPES, 100mM NaCl, pH 8.0.

Storage

Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Note

For research use only

Application notes

Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100 μg/ml. Dissolve the lyophilized protein in 1X PBS. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

Protein Sequence

MMASMRVVKE LEDLQKKPPP YLRNLSSDDA
NVLVWHALLL PDQPPYHLKA FNLRI SFPE
YPFKPPMIKF TTKIYHPNVD ENGQICLP II
SSENWK PCTK TCQVLEALNV LVNRPNIREP
LRMDLADLLT QNPELFRKNA EEFTLRFGVD
RPSLEHHHHH H

Purity

Greater than 95% as determined by reducing SDS-PAGE.

Source

E. coli

Uniprot ID

[O14933](#)

Entrez

[9246](#)

Expiration Date

6 months from date of receipt.