

Product Datasheet

Human UBE2R2 protein (orb392667)

Description

Recombinant human UBE2R2 protein is

and the target gene en...

Reactivity

Human

Endotoxins

Endotoxin: Less than 0.1 ng/μg (1 IEU/μg) as determined by LAL test.

Tested

ELISA, MS, SDS-PAGE, WB

Applications
Target

UBE2R2

Preservatives

Liquid in a 0.2 μM filtered solution of 50mM HEPES, 150mM NaCl, pH 7.5.

Form/Appearance

Liquid in a 0.2 μM filtered solution of 50mM HEPES, 150mM NaCl, pH 7.5.

Storage

Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Note

For research use only

Application notes

Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100 μg/ml. Dissolve the lyophilized protein in 1X PBS. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

Protein Sequence

MGSSHHHHHH SSGLVPRGSH MAQQQMTSSQ
KALMLELKSL QEPPVEGFRI TLVDESDLYN
WEVAIFGPPN TLYEGGYFKA HIKFPIDYPY
SPPTFRFLTK MWHPNYENG DVCISILHPP
VDDPQSGELP SERWNPTQNV RTILLSVISL
LNEPNTFSPA NVDASVMFRK WRDSKGKDKE
YAEIIRKQVS ATKAEAEKDG VKVPTTLAEY
CIKTKVPSND NSSDLLYDDL YDDDDIDDEDE
EEEDADCYDD DDSGNEES

Purity

Greater than 95% as determined by reducing SDS-PAGE.

Source

E. coli

Uniprot ID
Q712K3
Entrez
54926
Expiration Date

6 months from date of receipt.