

Product Datasheet

Human UBE2T protein (orb392669)

Description

Recombinant human UBE2T protein is

and the target gene enc...

Reactivity

Human

Endotoxins

Endotoxin: Less than 0.1 ng/μg (1 IEU/μg) as determined by LAL test.

Tested
Applications

ELISA, MS, SDS-PAGE, WB

Target

UBE2T

Preservatives

Liquid in a 0.2 μM filtered solution of 50mM HEPES, 150mM NaCl, 2mM DTT, 10% Glycerol, pH 7.5.

Form/Appearance

Liquid in a 0.2 μM filtered solution of 50mM HEPES, 150mM NaCl, 2mM DTT, 10% Glycerol, pH 7.5.

Storage

Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Note

For research use only

Application notes

Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100 μg/ml. Dissolve the lyophilized protein in 1X PBS. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

Protein Sequence

MGSSHHHHHH SSGLVPRGSH MQRASRLKRE
LHMLATEPPP GITCWQDKDQ MDDLRAQILG
GANTPYEKG VFKLEVIIPER YPFEPPIRF
LTPYHPNID SAGRICLDVL KLPPKGAWRP
SLNIATVLTS IQLLMSEPNP DDPLMADISS
EFKYNKPAFL KNARQWTEKH ARQKQKADDEE
EMLDNLPEAG DSRVHNSTQK RKASQLVGIE
KKFHPDV

Purity

Greater than 95% as determined by SEC-HPLC and reducing SDS-PAGE.

Source

E. coli

Uniprot ID

Q9NPD8

Entrez

29089

Expiration Date

6 months from date of receipt.