

Product Datasheet

Tubulin beta 1 Blocking Peptide (orb392744)

Description	The peptide is used to block Anti-Beta-tubulin antibody reactivity.
Reactivity	Human, Mouse, Rat
Tested Applications	Blocking
Target	TUBB1; TUBB2A; TUBB2B; TUBB4B; TUBB
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry. Experiment Notes: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q9BUF5 , P68371 , Q9CWF2 , Q13885 , P07437 , P68372 , Q99867 , Q7TMM9 , P85108 , Q13509 , Q9BVA1 , Q3ZCM7 , Q6P9T8 , Q3KRE8 , A2AQ07 , Q9H4B7 , P04350
Entrez	56604 , 498736 , 296554 , 7280 , 22151 , 291081 , 227613 , 84617 , 545486 , 10383 , 10382 , 10381 , 347733 , 203068 , 81027 , 347688 , 73710
Expiration Date	12 months from date of receipt.