

Product Datasheet

HDAC5 (Phospho-S259/220) Blocking Peptide (orb392773)

Description	The peptide is used to block Anti-Histone Deacetylase 5/9 (pS259/220) antibody reactivity.
Reactivity	Human, Mouse, Porcine, Rabbit
Tested Applications	Blocking
Target	HDAC5; HDAC9
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry. Experiment Notes: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q9UKV0 , Q99N13 , Q9UQL6 , Q9Z2V6
Entrez	79221 , 9734 , 10014
Expiration Date	12 months from date of receipt.