

Product Datasheet

KDELR2 Blocking Peptide (orb392792)

Description	The peptide is used to block Anti-KDELR2 antibody reactivity.
Reactivity	Human, Mouse, Porcine, Rat, Sheep, Zebrafish
Tested Applications	Blocking
Target	KDELR2
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry. Experiment Notes: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q9CQM2 , P33947 , Q5U305
Entrez	304290 , 66913 , 11014
Expiration Date	12 months from date of receipt.