

Product Datasheet

MAP3K5 (Phospho-S966) Blocking Peptide (orb392801)

Description	The peptide is used to block Anti-ASK1 (pS966) antibody reactivity.
Reactivity	Human, Mouse, Rat
Tested Applications	Blocking
Target	MAP3K5
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry. Experiment Notes: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
Purity	> 85%
Source	Synthetic
Uniprot ID	O35099 , Q99683
Entrez	26408 , 4217
Expiration Date	12 months from date of receipt.