

Product Datasheet

PAK1 (Phospho-S144/141/139) Blocking Peptide (orb392830)

Description	The peptide is used to block Anti-PAK1/2/3 (pS144/141/139) Antibody (orb393128) reactivity.
Reactivity	Human, Mouse, Rat
Tested Applications	Blocking
Target	PAK1; PAK2; PAK3
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q13177, Q61036, Q64303, Q13153, P35465, Q62829, Q8CIN4, O88643, O75914
Entrez	29433, 29431, 5063, 100910732, 18481, 5058, 5062, 224105
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)