

Product Datasheet

PRKCI (Phospho-T412/410) Blocking Peptide (orb392840)

Description	The peptide is used to block Anti-PKC iota/zeta (pT412/410) antibody reactivity.
Reactivity	Human, Mouse, Rabbit, Rat, Zebrafish
Tested Applications	Blocking
Target	PRKCI; PRKCZ
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry. Experiment Notes: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
Purity	> 85%
Source	Synthetic
Uniprot ID	F1M7Y5 , Q02956 , Q62074 , P41743 , Q05513 , P09217
Entrez	18759 , 18762 , 5584 , 84006 , 25522 , 5590
Expiration Date	12 months from date of receipt.