

Product Datasheet

PXN (Phospho-Y31) Blocking Peptide (orb392847)

Description The peptide is used to block Anti-Paxillin (pY31) antibody reactivity.

Reactivity Human, Porcine

Tested Applications Blocking

Target PXN

Form/Appearance Lyophilized powder

Storage Shipped at 4°C. Store at -20°C for one year.

Note For research use only

Application notes The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry. Experiment Notes: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.

Purity > 85%

Source Synthetic

Uniprot ID [Q66H76](#), [Q8VI36](#), [P49023](#)

Entrez [5829](#), [19303](#), [360820](#)

Expiration Date 12 months from date of receipt.