

Product Datasheet

alpha Synuclein (Phospho-Y133) Blocking Peptide (orb392869)

Description	The peptide is used to block Anti-Alpha-synuclein (pY133) antibody reactivity.
Reactivity	Human, Mouse, Porcine, Primate, Rabbit, Rat, Sheep
Tested Applications	Blocking
Target	SNCA
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry. Experiment Notes: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
Purity	> 85%
Source	Synthetic
Uniprot ID	O55042 , P37377 , P37840
Entrez	29219 , 20617 , 6622
Expiration Date	12 months from date of receipt.