

Product Datasheet

SP1 (Phospho-T739) Blocking Peptide (orb392870)

Description	The peptide is used to block Anti-SP1 (pT739) antibody reactivity.
Reactivity	Human, Primate, Rat
Tested Applications	Blocking
Target	SP1
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry. Experiment Notes: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
Purity	> 85%
Source	Synthetic
Uniprot ID	O89090 , P08047 , Q01714
Entrez	6667 , 20683 , 24790
Expiration Date	12 months from date of receipt.