

Product Datasheet

STMN1 (Phospho-S25) Blocking Peptide (orb392874)

Description	The peptide is used to block Anti-STMN1 (pS25) antibody reactivity.
Reactivity	Human, Mouse, Porcine, Primate, Rabbit, Rat, Sheep
Tested Applications	Blocking
Target	STMN1
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry. Experiment Notes: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
Purity	> 85%
Source	Synthetic
Uniprot ID	P13668 , P54227 , P16949
Entrez	29332 , 16765 , 3925
Expiration Date	12 months from date of receipt.