

Product Datasheet

TH (Phospho-S19) Blocking Peptide (orb392876)

Description	The peptide is used to block Anti-Tyrosine Hydroxylase (pS19) antibody reactivity.
Reactivity	Human, Mouse, Rat
Tested Applications	Blocking
Target	TH
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry. Experiment Notes: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
Purity	> 85%
Source	Synthetic
Uniprot ID	P04177 , P24529 , P07101
Entrez	7054 , 21823 , 25085
Expiration Date	12 months from date of receipt.