

Product Datasheet

ACO1 (Phospho-S711) Blocking Peptide (orb392898)

Description	The peptide is used to block Anti-Aconitase 1 (pS711) antibody reactivity.
Reactivity	Human, Mouse, Rabbit, Rat, Zebrafish
Tested Applications	Blocking
Target	ACO1
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry. Experiment Notes: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
Purity	> 85%
Source	Synthetic
Uniprot ID	P21399 , Q63270 , P28271
Entrez	11428 , 50655 , 48
Expiration Date	12 months from date of receipt.