

Product Datasheet

ATF2 (Phospho-T73) Blocking Peptide (orb392917)

Description	The peptide is used to block Anti-ATF2 (pT73) antibody reactivity.
Reactivity	Human, Mouse, Rat, Zebrafish
Tested Applications	Blocking
Target	ATF2
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry. Experiment Notes: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q00969 , P16951 , P15336
Entrez	11909 , 81647 , 102641666 , 1386
Expiration Date	12 months from date of receipt.