

Product Datasheet

CSNK1A1 (Phospho-Y294) Blocking Peptide (orb392961)

Description	The peptide is used to block Anti-CK1 alpha (pY294) antibody reactivity.
Reactivity	Human, Mouse, Rabbit, Rat, Sheep
Tested Applications	Blocking
Target	CSNK1A1
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry. Experiment Notes: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q8BK63 , P97633 , Q8N752 , P48729
Entrez	122011 , 113927 , 1452
Expiration Date	12 months from date of receipt.