

Product Datasheet

EPHA2 (Phospho-Y588/596) Blocking Peptide (orb392977)

Description	The peptide is used to block Anti-EPHA2 (pY588/596) antibody reactivity.
Reactivity	Human, Mouse, Primate, Rat
Tested Applications	Blocking
Target	EPHA2; EPHA3; EPHA4
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry. Experiment Notes: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
Purity	> 85%
Source	Synthetic
Uniprot ID	O08680 , P29320 , Q03145 , P29319 , P54764 , P29317 , Q03137
Entrez	13838 , 2042 , 1969 , 29210 , 2043 , 13836
Expiration Date	12 months from date of receipt.