

Product Datasheet

Estrogen Receptor beta (Phospho-S105) Blocking Peptide (orb392981)

Description	The peptide is used to block Anti-Estrogen Receptor beta (pS105) antibody reactivity.
Reactivity	Human, Mouse, Rat
Tested Applications	Blocking
Target	ESR2
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry. Experiment Notes: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
Purity	> 85%
Source	Synthetic
Uniprot ID	O08537 , Q92731 , Q62986
Entrez	13983 , 25149 , 2100
Expiration Date	12 months from date of receipt.