

Product Datasheet

GRIN2B (Phospho-Y1336) Blocking Peptide (orb393009)

Description	The peptide is used to block Anti-NMDAR2B (pY1336) antibody reactivity.
Reactivity	Human, Mouse, Rat
Tested Applications	Blocking
Target	GRIN2B
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry. Experiment Notes: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q00960 , Q01097 , Q13224
Entrez	24410 , 14812 , 2904
Expiration Date	12 months from date of receipt.