

Product Datasheet

HDAC7 (Phospho-S155) Blocking Peptide (orb393014)

Description	The peptide is used to block Anti-Histone Deacetylase 7 (pS155) antibody reactivity.
Reactivity	Human, Mouse, Rat, Zebrafish
Tested Applications	Blocking
Target	HDAC7
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry. Experiment Notes: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q8WUI4 , Q8C2B3 , Q99P96
Entrez	56233 , 51564
Expiration Date	12 months from date of receipt.