

Product Datasheet

HIST1H1B (Phospho-T17) Blocking Peptide (orb393016)

Description	The peptide is used to block Anti-Histone H1 (pT17) antibody reactivity.
Reactivity	Human, Mouse, Primate, Rabbit, Rat
Tested Applications	Blocking
Target	HIST1H1B
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry. Experiment Notes: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
Purity	> 85%
Source	Synthetic
Uniprot ID	P43274 , P16401 , P10412 , P43277 , P16402 , P43276
Entrez	14957 , 3007 , 50709 , 3008 , 56702 , 3009
Expiration Date	12 months from date of receipt.