

Product Datasheet

IL2RB (Phospho-Y364) Blocking Peptide (orb393033)

Description	The peptide is used to block Anti-CD122 (pY364) antibody reactivity.
Reactivity	Human, Mouse, Primate, Rat
Tested Applications	Blocking
Target	IL2RB
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry. Experiment Notes: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
Purity	> 85%
Source	Synthetic
Uniprot ID	P26896 , P14784 , P16297
Entrez	16185 , 25746 , 3560
Expiration Date	12 months from date of receipt.