


## Product Datasheet

### MST1R (Phospho-S1394) antibody (orb393114)

<b>Description</b>	Rabbit polyclonal antibody to MST1R (Phospho-S1394)
<b>Species/Host</b>	Rabbit
<b>Reactivity</b>	Human
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	WB
<b>Immunogen</b>	KLH-conjugated synthetic phosphopeptide corresponding to residues surrounding S1394 of human CD136 protein. The exact sequence is proprietary.
<b>Target</b>	MST1R
<b>Preservatives</b>	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
<b>Form/Appearance</b>	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
<b>Storage</b>	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
<b>Note</b>	For research use only
<b>Clonality</b>	Polyclonal
<b>Source</b>	Rabbit
<b>Uniprot ID</b>	<a href="#">Q04912</a>
<b>Entrez</b>	<a href="#">4486</a>
<b>Dilution Range</b>	WB: 1:500:1000
<b>Expiration Date</b>	12 months from date of receipt.



Western blot analysis of CD136 (Phospho-...