




Product Datasheet

MICA antibody (orb39313)

Description	Peptide affinity purified rabbit polyclonal antibody.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IF, IHC-P, WB
Immunogen	Synthetic Peptide
Target	MICA
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	42915
Uniprot ID	Q29983
NCBI	NP_000238.1
Dilution Range	IF/ICC: 1:25, FACS: 1:25, IHC-P: 1:25, WB: 1:1000
Expiration Date	12 months from date of receipt.



WB analysis of MDA-MB-231 cell line lysa...



WB analysis of (lane 1) A431 whole cell ...